Formula 1

Universal Investment for All Alloys and Pressable Ceramics

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

Formula 1 is a high performance phosphate investment designed to meet all your metal casting and glass pressing needs. It features both rapid and standard burnout as well as easy pouring and divesting capabilities. In addition, Formula 1 provides the widest expansion range and the best surface qualities of any investment available. To ensure success with Formula 1, please follow the instructions carefully.

Recommended Liquid: Special Liquid Concentrate – Plus

1. Pattern Preparation:

- Attach sprue to the wax pattern using pliable sticky wax.
- Mount patterns on crucible former base.
- Spray with Whip Mix SMOOTHX Wax Pattern Cleaner. Gently blow excess cleaner from pattern.

2. Preparation:

- If metal ring is used, line casting ring with a ceramic or dry paper liner positioned 6 mm short of each end. (Do not wet the ceramic liner.)
- Soak paper lined ring in water for a minimum of 1 minute. Then shake excess from ring before investing.
- Place ring onto former base.

3. Mixing:

For best results, store and use powder and liquid at room temperature (between 20° C (68° F) and 25° C (77° F)).

- Prepare liquid at suggested concentration following the chart below (distilled water recommended for dilution).
- Rinse bowl out with water and shake out excess. Use separate mixing bowls for phosphate and gypsum investments.
- Add measured liquid to mixing bowl. Incorporate powder by hand spatulation 10–15 seconds.
- Mechanical mix under vacuum, slow speed (350–600 RPM) for **2 minutes (120 seconds)**. Higher RPM mixers may require decreased mix time (90 seconds).

Note: Different paddle designs and mixer speed may require varying mixing time.

4. Benchset:

- Once the mold is poured, benchset 15 minutes minimum.
- Trim glaze off top of mold before burnout.
- For best results, place in preheated oven within 30 minutes of investing.

Note: Use rapid burnout method for benchset up to 1 hour. For longer than 1 hour, conventional burnout is recommended. Molds allowed to set more than 12 hours should be re-wet prior to burnout by soaking in water for 1–3 minutes.

5. Burnout:

Rapid Technique (Preheated Oven):

- Place molds in preheated oven at alloy manufacturer's recommended temperature, up to 925° C (1,700° F)*. For higher temperatures, place molds in oven at 925° C (1,700° F) then heat to final temperature at 14°-20° C (25°-35° F) per minute.
- Heat soak at final temperature for 30 minutes, add 10 minutes per additional mold.

***Note: Maximum preheat entry temperature for metal ring is 870° C (1,600° F).**

Standard Technique (Cold Oven):

- Place molds in oven at room temperature. Heat to desired temperature at 14°–20° C (25°–35° F) per minute.
- Heat soak at final temperature for 30 minutes, add 10 minutes per additional mold.

6. Casting (Alloy):

- Upon removal from oven, immediately cast according to alloy manufacturer's instructions.

Pressing (Pressable Ceramics):

- Press according to ceramic manufacturer's recommendations.
- Use 200 gram mold for restorations requiring two ingots.

7. Divesting:

- Allow metal castings and ceramic pressings to cool completely before divesting.

Note:

- Recommended Liquid: Special Liquid Concentrate – Plus
- For small volume mixes (less than 100 grams) decreasing liquid/powder ratio approximately 2 mL/100 gram will increase expansion and improve surface quality.
- To correct tight fits: Increase liquid concentration. If using metal rings, you may also use a double liner.
- To correct loose fits: Decrease liquid concentration. (Refer to Expansion Ratio Chart.)
- For large molds containing complex restorations or plastic sprues, runner bars or copings, the standard technique described above or a two-stage burnout technique is recommended.
- Special Liquid Concentrate – Plus may freeze during shipment in cold weather. If liquid should have crystallized upon arrival, it will no longer be usable. Keep liquid from freezing.

WARNING: Investments contain free silica — DO NOT BREATHE DUST. May cause delayed lung injury (silicosis/lung cancer).

Physical Properties:

Liquid/Powder Ratio	22 mL/100 gram
Working Time	6–8 minutes
Setting Expansion	1.6%
Thermal Expansion	0.65%
Compressive Strength	1,350 psi (9.2 MPa)

Expansion Ratio Chart

		TWO-MINUTE MIX TIME	LIQUID CONCENTRATION	60 GRAM 13 ML/60 GRAM		100 GRAM 22 ML/100 GRAM		150 GRAM 33 ML/150 GRAM		
				LIQUID (ML)	WATER (ML)	LIQUID (ML)	WATER (ML)	LIQUID (ML)	WATER (ML)	
EXPANSION	ALLOY	MORE	100%	13.0	0.0	22.0	0.0	33.0	0.0	
			90%	12.0	1.0	20.0	2.0	30.0	3.0	
		OPTIMUM	Base	80%	11.0	2.0	18.0	4.0	27.0	6.0
			Noble	75%	10.0	3.0	16.5	5.5	25.0	8.0
			High Noble	70%	9.0	4.0	15.0	7.0	23.0	10.0
	CERAMIC	LESS		60%	8.0	5.0	13.0	9.0	20.0	13.0
				50%	6.5	6.5	11.0	11.0	16.5	16.5
		MORE		90%	12.0	1.0	20.0	2.0	30.0	3.0
		OPTIMUM		80%	11.0	2.0	18.0	4.0	27.0	6.0
			Crowns, Veneers	70%	9.0	4.0	15.0	7.0	23.0	10.0
	Inlays, MODs	60%	8.0	5.0	13.0	9.0	20.0	13.0		
	LESS		50%	6.5	6.5	11.0	11.0	16.5	16.5	
			40%	5.0	8.0	9.0	13.0	13.0	20.0	

*Recommended concentrations are approximate and can be adjusted to optimize fit. To increase expansion, use more liquid and less water. To decrease expansion, use less liquid and more water. Always maintain total liquid / water volume.

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