



DIGITAL WAX CARVING PENCIL – DUAL

OPERATION MANUAL



Whip Mix Corporation

P.O. Box 17183 • Louisville, KY USA 40217-0183

502-637-1451 • 800-626-5651 • Fax 502-634-4512 • www.whipmix.com

IMPORTANT SAFEGUARDS



**CAUTION: FAILURE TO FOLLOW THESE SAFEGUARDS
MAY RESULT IN FIRE OR ELECTRIC SHOCK.**

1. Make sure that the unit is connected to a power source that is grounded.
2. To prevent burning, do not touch unit when in use.
3. Never touch unit, pencil, cord or plug with wet hands.
4. Protect the unit from bumps or scrapes, which may cause mechanical problems.
5. Do not use the unit if the power cord or plug is damaged in any way.
6. Make sure that the power cord is not twisted, bent or pinched.
7. To prevent possible overheating, do not block the ventilation openings at the bottom of the unit.
8. Unplug the unit during storms and when it will not be used for an extended period.
9. Do not pull on the cord when unplugging the unit.
10. Do not store or use the unit in a room with high humidity or near dust or spraying water.
11. Keep unit away from heating instruments and other heat sources.
12. Always unplug the unit before cleaning. Do not spray water directly on the unit. Avoid chemical cleaning products.
13. Keep power plug prongs free of dust and water.
14. If you detect noise or smell smoke coming from the unit, unplug immediately.
15. Do not disassemble or attempt to repair the unit.
16. Turn the unit off and let cool before changing carver tip.
17. Do not allow tip of carver to touch the wire.

Please contact Whip Mix with any questions or for troubleshooting advice.

Use of the equipment in any way not intended by the manufacturer compromises the safety of the equipment.

OPERATING INSTRUCTIONS

1. Plug in unit and press the ON/OFF switch. Green Power switch will illuminate.
2. Insert desired tip into pencil and insert cord into “OUT 1” or “OUT 2” jack on unit.
3. Set the desired temperature using the UP/DOWN button. Press the button once to raise or lower the temperature by one degree or hold the button down for quicker temperature setting.
4. Press C/F to select either Celsius or Fahrenheit scale. (C/F will appear in display window.)
5. If an error occurs, a buzzer will sound and the window will display ERR. Turn Power switch off and then back on to reset.

FUNCTIONS & FEATURES

- Easy to use.
- Excellent heat conductivity.
- Long service life.
- Easy-to-use, visible control buttons.
- Two-line circuit, designed for simultaneous work.
- Minimal temperature deviation.
- Accurate temperature maintenance with reset function.
- Built-in digital CPU circuit accurately maintains desired temperature.
- Digitally displays temperature in either Celsius or Fahrenheit.
- Previous working temperature is kept in memory, even if the unit has been turned off.
- Usable with AC100-240V without additional adjustment .
- Silicone handle reduces hand fatigue.
- Convenient magnetic cradle on surface or external cradle.
- Translucent compartment in rear of unit for easy storage of carving tips.

TECHNICAL SPECIFICATIONS

Voltage: AC 100-240V 50/60Hz 20W

Dimensions: 6.0"W × 4.9"D × 2.9"H
152 mm × 123 mm × 70 mm

Temperature:

Optimum Indoor Temperature: 5°C – 40°C (41°F – 104°F)

Optimum Indoor Humidity: 50% – 80%

Weight: 2.8 lbs. (1.27 kg)



**NOTE: MAIN VOLTAGE SUPPLY SHOULD NOT
EXCEED ± 10% OF THE NOMINAL VALUE**

PACKING LIST

1. Control Box
2. 2 Electric Pencils
3. 2 Pencil Stands
4. AC Cord
5. Two (2) Tips

Description	Item #
Digital Wax Carving Pencil–Dual	10150



more reasons. one source

P.O. Box 17183 • Louisville, KY USA 40217-0183
502-637-1451 • 800-626-5651 • Fax 502-634-4512
www.whipmix.com