ProCal™ User Manual

Product Overview

The Pro Cal Calibration Kit is a professional grade instrumentation system designed to rapidly monitor the calibration status of your furnace, and automatically recalibrate the furnace to factory specifications.

It is recommended that you read this manual thoroughly before using the kit as there are valuable examples of procedures and decisions you will want to be acquainted with in order to maximize the benefits you've invested in.

What You Can Do With It

QC = Quick Check (Pro 200 Series). For FireLite refer to page 9.

The Quick Check feature allows you to rapidly spot check the need for calibration. Simply place the Special Firing Tray on the Lift Table, connect the Tray to the Instrument Box and the Instrument Box to the furnace calibration port, and dial in the desired temperature you want to verify. The LCD readout on the ProCal box should agree with the furnace display to within three degrees Celsius. If not, you can choose one of several options:

1. If you check several temperatures and only one seems to be different, you can go into the firing programs which **have hold temperatures at that value** and use the **Program Level Fusing Adjust** to fine tune just those programs.
2. If you check several temperatures and **they all either overfire or underfire**, you can go into the **Special Functions** menu and select **High Fusing Adjust for temperatures over 800°C** or **Low Fusing Adjust for temperatures under 800°C** and enter the average correction value needed.
3. You can do a fully automated calibration.

Full Recalibration (Pro 200 Series Only)

The full recalibration mode requires several hours to complete, so you'll want to reserve it for times when the temperature is **more than three degrees** from nominal. Why does it take so long? To get a valid calibration which will most closely match the conditions in which it is used every day the furnace must be fully heat soaked. Secondly, the automatic calibration program requires one pass at three separate temperatures to build the calibration file. Thirdly, a second complete pass is necessary to confirm that the calibration pass is valid. This is exactly the same rigorous process that your furnace went through before it left the factory originally.

Why Do Furnaces Require Periodic Temperature Recalibration?

There are several reasons why furnaces require periodic QC checking and recalibration:

Effluents given off by the materials fired in the furnace can deposit on the platinum/palladium thermocouple and

in some cases, chemically react with the noble metals, thereby affecting the temperature sensed.

Oxide layers can build up on the thermocouple through many cycles of vacuum, then exposure to air and chemically react with the noble metals in the thermocouple.

Prolonged high heat exposure can eventually crystallize the noble metals making them brittle and inhibiting the buildup of charge on each leg. The millivolts given off by a thermocouple to register temperature accumulate all along the heat exposed legs of the thermocouple, not just at the tip.

The heater wire also becomes brittle and scaly along its surface after extended usage. This can affect the infrared spectrum irradiated by the heater and redirect the angle which the photons emitted by the heater such that the ones striking the thermocouple are more diffused and therefore create less charge buildup in the thermocouple, affecting the calibration.

How It Works - What's it doing in there?

The thermocouple in the Special Firing Tray monitors the furnace activity and reports its accumulated charge back to the Instrument Box. The Instrument Box filters the electrical noise inherent with very low level signals in the microvolt to millivolt range. The cleaned-up signal is then sent to a pre-amplifier with a built in digital filter to further boost the signal and clean it up some more. Finally, the signal is presented to a 24 bit Analog to Digital Converter (ADC) for conversion into digital ones and zeros to be transmitted to the furnace computer. There is also a local "ambient" temperature sensor inside the Instrument Box which monitors the "cold junction" temperature at the point where the thermocouple alloys transition into the electronics. This cold junction temperature sensor output is automatically used by **the computer inside the Instrument Box** to correct any measurement errors as the reference signal from the Special Firing Tray makes its way back to the furnace.

Once the digital reference value from the Special Firing Tray is received back in the furnace computer its value is compared to the internal furnace thermocouple and the sign (relative over-firing or under-firing) and the magnitude of the difference are stored permanently in the furnace computer.

After the calibration process is complete, every time the furnace temperature is commanded to a specific temperature, the furnace computer looks at the target temperature, determines which calibration value is appropriate for that target value, and either adds or subtracts the stored calibration correction value in order to produce the accurate temperature desired, all in real-time.

Kit Contents

- Instrument Box
 - Special Firing Tray
 - Data Cable
 - Operator Manual
 - FireLite 95031 Wall Charger (Sold separately)
-

Important Notes / Reminders

The Special Firing Tray and the Instrument box form a **MATCHED SET**. These two components are calibrated together at the factory and the overall system accuracy depends upon these two components always being used together. If more than one Kit is on site, **never use the Special Firing Tray from one Kit with the Instrument box of another Kit or inaccurate calibrations will result.**

Since High Temperatures are present in the Furnace for extended periods of time, **it is NOT RECOMMENDED that Full Calibrations be attempted unless someone is present in the vicinity of the furnace at all times.**

Component Familiarization

Special Firing Tray

Instrument Box

Data Cable

FireLite 95031 Wall Charger (optional)


Setup for Either Quick Check or Full Calibration (Pro 200 only)

Turn the Furnace Power to the OFF position.

This is required to prevent electrical surges which could damage the sensitive electronics in the Instrument Box.

Remove the two small thumbscrews on the right side of the furnace chassis to expose the calibration ports. There are two 9-pin connectors inside the furnace which are now exposed.


Turn the Furnace Power to the ON position.

Wait for the furnace self-test to complete and for the Lift Table to automatically lower to the bottom position.

• Tray Placement

Carefully remove the normal furnace firing tray by lifting it about an inch upward to clear the centering pin with a pair of forceps or tongs.


Use the forceps or tongs to place the Special Firing Tray with the thermocouple embedded in it onto the tray centering pin and allow it to come to rest on the Lift Table. Swivel the Special Firing Tray so that the wire umbilical coming out of the edge of the Special Firing Tray is facing the front of the furnace.

Use the Furnace **Up** Arrow button to raise the table until it is fully closed and the wire umbilical is protruding out of the front/center of the muffle.

NOTE: IT IS IMPOSSIBLE TO PULL VACUUM WITH THE SPECIAL FIRING TRAY IN PLACE BECAUSE OF THE WIRE UMBILICAL. THE CALIBRATION QUICK CHECKS AND FULL RECALIBRATION DO NOT REQUIRE VACUUM, SO THIS IS NOT AN ISSUE.

- Finally, connect the end of the Special Firing Tray umbilical into the Instrument box by inserting the green connector into the end of the Box opposite the end with the Data Cable plugged into it.

IMPORTANT: PLACE THE INSTRUMENT BOX WELL TO THE RIGHT SIDE OF THE FURNACE SUCH THAT THE BOX IS OUT OF THE PENUMBRA OR "SUNSHINE" EMITTED BY THE MUFFLE HEATER. IT NEEDS TO BE "IN THE SHADE" SO THAT THE ELECTRONICS ARE NOT AFFECTED BY DIRECT HEAT "SHINING" ON THE INSTRUMENT BOX. (EVEN IF YOU HAVE TO LET THE INSTRUMENT BOX HANG OFF THE FRONT OF THE TABLE OR COUNTER SUPPORTING THE FURNACE.)

Quick Check Survey (Assumes the Physical Setup Step is Complete)

Procedure for Pro 200 Series Only

Navigating to the Pro Cal Menu

From the Pro 200 Main Menu Select **Special Functions**.

Main Menu
>Special Functions

Scroll Through the **Special Functions** Menu Items until you reach “Calibrate Oven” and select it.

>Calibrate Oven

When prompted for a password , enter “291”

Password = 291

which will Bring up the **Pro Cal** Menu as below.

Pro Cal Furnace
>Quick QC Check
Pro Cal - Cool Muffle
Pro Cal - Soaked Muffle

Select “Quick QC Check” by pressing the **Selector Knob** OR pressing the Next button.

Entering a Temperature to Check

The Pro 200 is now ready to raise the temperature to any setting you desire and hold the temperature constant for up to 40 minutes. Select the specific temperature by rotating the **Selector Knob** clockwise for increased temperature or counterclockwise for decreased temperature. Once rotated to the desired temperature press the **Selector Knob** inward to got to the newly selected temperature.

Pro Cal Furnace
Change Idle Temp
>Roll Knob 250 C
Press Knob = GoTemp

NOTE: YOU CAN SELECT A NEW STABLE IDLE TEMPERATURE AT ANY TIME BY PRESSING THE LAST BUTTON ON THE PRO 200 ROLLING THE SELECTOR KNOB TO A NEW TEMPERATURE VALUE AND THEN PRESSING THE SELECTOR KNOB AGAIN. THE PRO 200 WILL TRANSITION TO THE NEW VALUE. YOU CAN REPEAT THIS AS MANY TIMES AS YOU WISH. IF YOU ARE DONE WITH THE PRO CAL QUICK CHECK ACTIVITY YOU CAN PRESS THE MAIN MENU BUTTON TO EXIT THE PRO CAL MODE AT ANY TIME AND THE PRO 200 WILL AUTOMATICALLY REVERT BACK TO THE PRIMARY IDLE TEMPERATURE YOU DESIGNATED WHEN YOU ORIGINALLY SETUP YOUR FURNACE FOR THE FIRST TIME.

The Pro Cal “Idle” temperature is NOT the same as the Idle Temperature you selected from the Special Functions

Menu. The Pro Cal is a separate operating mode.

Allow Stabilization Time

It is suggested that a period of at least 15 – 20 minutes is allowed for stabilization time before manually recording the value on the Instrument box LCD display to insure a stable reading.

What to do with the Results

It is suggested that you survey three temperatures (minimum) to look for overall trends, and to help you decide which type of correction method best fits your available time and addresses the current situation the best.

The three primary temperatures suggested are:

Celsius Value	Fahrenheit Value
850°C	1562
1000°C	1832
1150°C	2102

These are the temperatures used by the Full Calibration method.

For example, if you surveyed the three temperatures above and when the furnace display said 850°C the LCD display on the Instrument box read 851°C.

AND at 1000°C on the furnace display the Instrument box read 1002°C.

AND at 1150°C on the furnace display the Instrument box read 1151°C.

All three being slightly higher on the Instrument box indicates that the furnace is slightly OVERFIRING across the board. Since the Over-firing values are not more than 3 degrees C above the nominal target temperatures, a Full Calibration run IS NOT INDICATED. The fastest and easiest way to proceed from here is to enter a 1 or 2 degree OVERFIRE into the **Special Functions** High Fusing Adjust section of the Pro 200 software.

The procedure for doing so follows in the next section of this manual.

One thing to remember about using the **Special Functions** Menu selected **Fusing Adjustment** entry is it affects all User Program Firings equally, every single one of them.

Conversely, if you had found, in the QC Survey above:

Furnace Display	Instrument Box Display
850°C	850°C
1000°C	1001°C
1150°C	1150°C

Only a single temperature range is over-firing. You would NOT want to use the **Special Functions** High Fusing Adjust entry method to register the Over-firing at the

1000°C temperature because the following would result:

Furnace Display	Instrument Box Display
850°C	849°C
1000°C	1000°C
1150°C	1149°C


As you can see from the table above, using the High Fusing Adjustment in this instance would fix the 1000°C Overfire concern, but it would also produce an undesired Under-firing situation in the ranges above

(1150°C) and below (850°C).

You can still address a deviation from nominal in a single range by going into each program that has a **Hold Temperature** in that range and use the **Program Level Fusing Adjustment**.

How To Enter An Overfire Correction Using The Special Functions High Fusing Adjust (Pro 200 only)

(Affects all User Firing Programs with a Hold Temperature ABOVE 800°C)

Navigate to screen #3 above. If you want to register an Over-firing correction, simply press the **Selector Knob** OR press the **Next** button. If you want to register an Under-firing correction, first roll the **Selector Knob** until the answer to the question “High Overfire?” displays “No” (meaning its not Over-firing, its Under-firing), then press the **Selector Knob** OR press the **Next** button to get to screen #4 above.


Enter the difference between what the Furnace Display said during the Quick Check and what the Instrument box Display said at the same time. Then press the **Selector Knob** OR press the **Next** button to record the correction.

Example: Lets say the furnace Target Temperature was set to 1000°C. When stabilized, the furnace Displayed 1000°C, and the Instrument Box Display shows 998°C. **The Instrument Box is always the reference point, not the furnace.** So, the furnace is running two degrees **HIGHER** than the Instrument Box, in other words, the furnace is **OVERFIRING** by two degrees. You would answer “Yes” to the “High Overfire?” question in screen #3 above (Yes is the default already displayed) by pressing the **Selector Knob** OR pressing the **Next** button, and then in screen #4 above you would roll the **Selector Knob** until the adjustment value “2C” displays, then press the **Selector Knob** OR press the **Next** button to record the correction.

Remember, this correction will affect all programs with Hold Temperatures above 800°C.

How to Enter An Overfire Correction Using The Program Level Fusing Adjust (Pro 200 only)

Affects This One User Firing Program ONLY)

Full Calibration Run (Pro 200 only)

IMPORTANT NOTE: WHEN THE FULL CALIBRATION RUN FINISHES SUCCESSFULLY, THE PRO CAL SOFTWARE WILL AUTOMATICALLY ZERO OUT ANY PREVIOUSLY ENTERED HIGH FUSING ADJUSTMENTS AND LOW FUSING ADJUSTMENTS AS THE NEWLY MINTED FULL CALIBRATION CANCELS THE NEED FOR SUCH ENTRIES. THE PROGRAM LEVEL FUSING EDITS WILL, HOWEVER REMAIN IN EFFECT UNTIL USER MODIFIES THEM AGAIN.

When to Run

It is recommended that a Full Calibration Run be performed **ONLY** when the QC Quick Check survey indicates a deviation of over 3 degrees C between the displayed value on the furnace display and the displayed value on the Instrument Box display.

To Preheat or Not

It is recommended that the Preheat version of the Full Calibration Run be used for best accuracy. The difference in Run Time is 27% but even the shorter version requires 2:15 hours versus the cool muffle version Run Time of 4:25 hours.

How to Launch the Calibration Run

See the Physical Setup Section Described Previously to configure the Kit for a Run, then:

From the Pro 200 **Main Menu** Select **Special Functions**.

Main Menu

>Special Functions

Scroll Through the **Special Functions** Menu Items until you reach “Calibrate Furnace” and select it.

>Calibrate Furnace

When prompted for a password , enter “291” which will Bring up the **Pro Cal** Menu as below.

Password = 291


Scroll down to the Pro Cal with Cool Muffle Selection as shown below and press the **Selector Knob**.


Pro Cal Furnace
Quick QC Check
>Pro Cal – Cool Muffle
Pro Cal – Soaked Muffle

How to Keep Tabs on the Calibration Run


The sample displays below will keep you informed of the run’s progress

COOL MUFFLE TIMELINE

SOAKED MUFFLE TIMELINE**Inspecting the Results**

- **Within 40 Minutes of Completion**
The three calibration results for 850°C, 1000°C, and 1150°C are on the confirmation screen shown above.
- **Anytime Thereafter** you can retrieve the information by doing the following:

Detailed Directions for Use with the FireLite Oven

1. Replace the firing tray with the Pro Cal Tray.
2. Power on the FireLite oven and set the idle temperature (Special Function 22) to 650 degrees C/1202 degrees F, and let it soak for at least an hour with the table down.
3. Plug the connector from the Pro Cal Tray into the Pro Cal device.
4. Plug the FireLite power adaptor into the Pro Cal device and into an AC wall power outlet (115 Volt ONLY). Make sure the Pro Cal device display powers on.
5. Determine what temperature needs to be checked and create a program with that temperature as the Final Temperature with a Hold Time of 30 minutes.
Note: For High Fusing porcelains select a temperature ABOVE 800 degrees C.
For Low Fusing porcelains select a temperature BELOW 800 degrees C.
6. Run the program and 15-20 minutes into the Hold Time record the difference between the oven temperature and the Pro Cal device displayed temperature.
7. Abort the program and enter the Fusing Adjust mode:
(Special Functions 18 for HIGH FUSING OVER-FIRES)
(Special Functions 19 for HIGH FUSING UNDER-FIRES)
(Special Functions 16 for LOW FUSING OVER-FIRES)
(Special Functions 17 for LOW FUSING UNDER-FIRES)
8. Enter the recorded value for the difference found in Step 6 during the program run, press Enter to save.
9. This completes the Fusing Adjust process.

Troubleshooting (Pro 200 only)

- Possible Error Conditions & What to Do

Error Number Displayed	Problem	What to Do About It
100	Pro Cal reading exceeds upper limit	This could be a bad furnace thermocouple, causing the furnace temperature to go very high. Remove the Pro Cal Kit and verify the furnace Actual Temperature is tracking the furnace Target Temperature. If the furnace itself is not the issue, contact Whip Mix Customer Service.
102	The difference between the furnace temperature reading and the Pro Cal temperature reading is too great.	Make sure the Pro Cal tray is inside the furnace, and check the cable connections between the tray and the Instrument Box and the Instrument Box to the furnace connectors are seated tightly. Restart the Pro Cal cycle, if error persists, contact Whip Mix Customer Service.
103	Cannot read temperature from Instrument Box.	Make sure the Pro Cal tray is inside the furnace, and check the cable connections between the tray and the Instrument Box and the Instrument Box to the furnace connectors are seated tightly. Turn the furnace main power off, leave it off for ten seconds then turn it back on and restart the Pro Cal cycle, if error persists, contact Whip Mix Customer Service.
104	Cannot sense temperature rising on Pro Cal Tray (but electronic communications are good)	This most likely means the Pro Cal Kit is connected correctly, but the Pro Cal tray is not physically on the furnace lift table. Verify Pro Cal Tray is inside furnace muffle. If error persists, contact Whip Mix Customer Service.
105	High Ambient Temperature Reading	The Pro Cal Instrument Box contains an "ambient temperature sensor" which monitors the internal temperature of the Instrument Box itself. If the reading is too high it can affect the accuracy of the Calibration. Move the Instrument Box away from the furnace so that it is not exposed to elevated heat from either the furnace being calibrated or another nearby furnace.

Troubleshooting

(Pro 200 only)

106	Unstable temperature readings.	The primary cause would be an intermittent electrical connection between the Tray and the Instrument Box or between the Instrument Box and the furnace. Re-seat all connections and re-launch the calibration. Turn the furnace main power off, leave it off for ten seconds then turn it back on and restart the Pro Cal cycle, if error persists, contact Whip Mix Customer Service.
107	No Ambient Temperature Reading	(See also Error 105 above) This is most probably an electronic fault inside the Instrument Box and would render inaccurate calibrations if not addressed. Turn the furnace main power off, leave it off for ten seconds then turn it back on and restart the Pro Cal cycle, if error persists, contact Whip Mix Customer Service.

The following three Error Codes 108, 109 and 110 have multiple possibilities.

First, verify all the Pro Cal Kit electrical connections are tightly connected, especially if Error 108 is displayed (the first temperature to be checked). Re-start the Pro Cal cycle afterwards.

Second, if the furnace has been idle for a while or at low heat with the Lift Table open, moisture may have accumulated in the white muffle insulation. Pre-soak the muffle at high heat for 30 minutes to an hour (you can use the Quick Check Selection on the Pro Cal Menu and enter 850 – 1000°C as a target temperature), then re-start the Pro Cal cycle afterwards.

Third, use the Quick Check mode and set an 850°C Target Temperature and wait for the furnace Actual Temperature to rise to 850°C. Observe the Actual Temperature for a couple of minutes and if the Actual Temperature on the furnace display is not stable at 850°C, there is a significant problem in the furnace itself which will require attention. Consult Whip Mix Customer Service in this instance.

Error Number Displayed	Problem	What to Do About It
108	Pro Cal Temperature Check failed at 850°C	See above.
109	Pro Cal Temperature Check failed at 1000°C.	See above
110	Pro Cal Temperature Check failed at 1150°C.	See above

These are two additional Error Codes.

Error Number Displayed	Problem	What to Do About It
111	Communication Error	The furnace cannot communicate with the Pro Cal Instrument Box electronically. Check the connections at both the furnace and the Instrument Box. (Possible defective cable.) Turn the furnace main power off, leave it off for ten seconds then turn it back on and restart the Pro Cal cycle, if error persists, contact Whip Mix Customer Service.
112	Unreliable temperature data readings.	Basically the same as Error 106. Turn the furnace main power off, leave it off for ten seconds then turn it back on and restart the Pro Cal cycle, if error persists, contact Whip Mix Customer Service.

Technical Support


For additional support please contact us at (800)626-5651.

Specifications

- Special Firing Tray – Type K Thermocouple with ceramic support, ceramic fiber sheath in hot zone.
Maximum Temperature: 1200 Degrees Celcius.
- Instrument Box – ABS plastic enclosure with 2 line by 8 character LCD Display.
Maximum Temperature: 125 Degrees F
Power: +5 VDC +/- .1% @ .030 Amps, furnished by Furnace under test via Data Cable.
Fuse: None (draws power from furnace under test)
- Data Cable – 9 pin DSub Male/ 9 Pin DSub Female, shielded. Standard RS232 Computer Cable.
Maximum Temperature: 125 Degrees F.
- System Accuracy: +/- 0.125% of measured temperature.
- Calibration Results Accuracy: +/- 3 degrees Celsius from 700°C to 1200°C.
- Pollution Category: 1
- For Indoor Use Only
- Environmental Ratings:
Altitude: Up to 2,000 meters
Ambient Temperature: 5°C – 40°C
Relative Humidity: 80% @ 31°C
linearly decreasing to 50% @ 40°C
- Duty Cycle: 100%
- Dimensions: Case – 14.0", (355 mm) Width
14.5", (369 mm) Length
5.0", (127 mm) Depth
- Shipping Weight: 3.86 LB, (1.75 Kg)

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