



DENTAL TECHNOLOGY SOLUTIONS

# Ongoing Whip Mix Print Resin / 3D Printer Qualification















Qualification Complete = ▲ Ongoing Qualification Process = ●	VeriModel OS Ivory (385/405nm)	VeriModel OS White (385/405nm)	VeriModel OS Grey (385/405nm)	VeriModel OS Golden Brown (385/405nm)	VeriTray (385nm)	VeriCast (385nm)	VeriGuide OS Clear (385/405nm)	VeriSplint OS (385nm) (510K cleared)	Dentca	Surgical Guide for Carbon
 Whip Mix VeriBuild	▲	▲	▲	▲			▲		▲	
 Asiga MAX (385nm/405nm)	▲	▲	▲	▲	▲	▲	▲	▲	▲	
 Asiga Pro2 75 (385nm/405nm)	▲	▲	▲	▲	▲	▲	▲	▲	▲	
 Ackuretta Freeshape 120 (405nm)	▲	▲	▲	▲			▲			
 Ackuretta Diplo (405nm)		▲	●	▲			▲			
 Anycubic (405nm)		▲		▲			▲			
 FlashForge Hunter (405nm)		▲		▲			▲			
 Microlay Versus (385nm/405nm)	●	▲		●	●	●	▲	●		
 NewPro3D NP1 (385nm/405nm)	●	▲		▲	▲	▲	●	●		

Continued on reverse side



DENTAL TECHNOLOGY SOLUTIONS

# Ongoing Whip Mix Print Resin / 3D Printer Qualification

Qualification Complete = ▲ Ongoing Qualification Process = ●	VeriModel OS Ivory (385nm) 	VeriModel OS White (385/405nm) 	VeriModel OS Grey (385/405nm) 	VeriModel OS Golden Brown (385/405nm) 	VeriTray (385nm) 	VeriCast (385nm) 	VeriGuide OS Clear (385/405nm) 	VeriSplint OS (385nm) (510K cleared) 	Dentca 	Surgical Guide for Carbon 
 <b>Carbon 3D M1, M2, M2D</b>									▲	▲
 <b>Makex M-dental</b>		▲					▲	●		
 <b>Makex M-one pro</b>					▲	▲				
 <b>Photocentric LC Dental</b>					▲	▲	▲	●		

Don't see your printer? Contact us at [www.whipmix.com/contact-us](http://www.whipmix.com/contact-us) Visit [www.whipmix.com](http://www.whipmix.com) for updates.

Whip Mix Corporation • 361 Farmington Avenue • Louisville, KY USA 40209  
 Toll-Free: 800-626-5651 • Phone: 502-637-1451 • [www.whipmix.com](http://www.whipmix.com)