



# SCREENING SYSTEM COMPARISON

	Nonin Pulse Oximeter / FoxPro Screening System	Nox T3 HST Screening System
		
Number of Channels	2	12
<b>Data Captured:</b>		
Oxygen Saturation (per hour)	■	■
Oxygen Desaturation Index	■	■
Heart Rate Range	■	■
Apnea - Hypopnea Index		■
Audio (Snoring and Tooth Grinding)		■
Snoring Index		■
Body Position		■
Bruxism Episode Index		■
Break of Bruxism Episodes (Clench, Grind)		■
<b>Hardware Included:</b>		
Monitor		■
Pulse Oximetry	■	■
Software	■	■
Consumables	Batteries	■
Charging Cables	■	■
Case		■
<b>Reporting Capabilities:</b>		
Wellness Report	■	
Dental Bruxism Report		■
Physician Interpretation		■

**Want more information? Contact us at [www.whipmix.com/contact-us](http://www.whipmix.com/contact-us)**

Whip Mix Corporation ■ 361 Farmington Avenue ■ Louisville, KY USA 40209

Toll-Free: 800-626-5651 ■ Phone: 502-637-1451

[www.whipmix.com](http://www.whipmix.com)