



Go Green With Envy



Easy on you

Easy on instruments

Easy on the environment

Many cleaning products contain active ingredients, such as sodium carbonate, which may cause burns and irritate the eyes. Prolonged immersion of metal trays in these solutions may also cause corrosion. Furthermore, release of these hazardous chemicals into the water system may be toxic to aquatic plants and animals.

Cleans without harsh chemicals

Green Envy™ dissolves alginate and plaster residues by interaction rather than by attacking them with corrosive chemicals. The residues rinse away—no need to scrape or scour.

pH neutral

Green Envy™ is neither acidic or alkaline and is gentle to the skin. It can be safely handled and will not cause corrosion to metal instruments.

Non-toxic

Green Envy™ does not contain toxic compounds that pose a risk to users. This also reduces hazards in the municipal waste stream and damage to the ecosystem from persistent chemicals.

Bio-degradable

Green Envy™ is broken down by micro-organisms into simpler, more stable compounds. It decomposes without causing environmental pollution.

Easy to use

Simply mix Green Envy™ with water according to the directions, immerse the instruments for 15-30 minutes, then brush or rinse with warm water. The only thing left behind is a fresh, light citrus scent.

Green Envy™ **TRAY CLEANER**

Item No. 09635
2.2 lb. (1 kg) Jar

Call Whip Mix at
800-626-5651
or your dental dealer.



RESTORATIVE ORAL HEALTH
Whip Mix Corporation
361 Farmington Avenue
Louisville, KY 40209
www.whipmix.com