

### Ongoing Whip Mix 3D Print Resins /Post Curing Units Qualifications

3D print resin Curing unit	Verimodel OS White	Verimodel OS Ivory	Verimodel OS Golden Brown	Verimodel OS Grey	VeriGuide OS	VeriSplint OS	VeriCast	VeriTray	Dentca
<b>G171/Otoflash</b>	2X 2000 flashes	2X 2000 flashes	2X 2000 flashes	2X 2000 flashes	6000 flashes	6000 flashes	1 x 750 flashes	2X 2000 flashes	Pending
<b>IntelliRay 600W/Uvitron</b>	pending	pending	pending	pending	5 min at the lowest settings of the rack	5 min at the lowest settings of the rack	pending	pending	10 minutes per side @ 50% intensity in Glycerin bath
<b>SunRay 400W /Uvitron</b>	pending	pending	pending	pending	5 min at the lowest settings of the rack	5 min at the lowest settings of the rack	pending	pending	10 minutes per side in Glycerin bath