



# Cadiax<sup>®</sup> Compact 2

## Mandibular Recording Device

The Cadiax<sup>®</sup> Compact 2 System is designed for registration and display of hinge axis movements of the human mandible. The data collected is used to program the articulator providing more predictable restorative results.

- ▶ Patient set up and recordings are accomplished in less than 10 minutes.
- ▶ Easily incorporated into any teaching or practice philosophies.
- ▶ Visual aid and screening tool for TMJ anomalies.



RESTORATIVE ORAL HEALTH

# Cadiax<sup>®</sup> Compact 2

## Mandibular Recording Device

### New Features

- ▶ Real time recording to visually locate the reference position before recording.
- ▶ Increased resolution.
- ▶ Compatible with Vista, Windows XP, and Windows 2000.
- ▶ Optimized for recording biomechanical movements of the TMJ through new magnetic technology.
- ▶ Multiple recordings offer:
  - Comparison of tracings
  - Verify accuracy
- ▶ New streamlined Cadiax unit:
  - Space saving dimensions  
1.5" w x 6.5" h x 6.5" d
  - Portable
  - Plugs directly into computer  
USB port

Developed and produced by Gamma Med. Wiss. Fortb. GmbH, Austria. [www.gammadental.com](http://www.gammadental.com)

Distributed by:  
Whip Mix Corporation - West  
1730 E. Prospect Road, Suite 101  
Fort Collins, CO 80525  
800-626-5651  
[www.whipmix.com](http://www.whipmix.com)

### Cadiax<sup>®</sup> Compact 2 System and Recording Software

Item 20002375



### Gamma Dental Software<sup>®</sup>

Item 20002916



RESTORATIVE ORAL HEALTH