

INNOVATIVE TECHNOLOGY



DIGITAL

---

# Waxing Units

---

HELPING SHAPE BETTER DENTAL RESULTS



*more reasons. one source*

## DIGITAL WAX CARVING PENCILS

### Standard Features

- Easy to use
- Excellent heat conductivity
- Long service life
- Accurately maintains desired temperature
- Digitally displays temperature in either Celsius or Fahrenheit temperature scales
- Previous working temperature is stored in memory, even if the unit has been turned off
- Silicone cool-grip handle reduces hand fatigue and temperature transfer
- 2-year warranty



### Digital Wax Carving Touch Pencil

#### In addition to standard features:

- Touch-controlled temperature setting
- Two-line circuit, designed for simultaneous work
- Convenient magnetic cradle on surface
- Compartment in rear of unit for easy storage of carving tips
- Comes complete with six of the most commonly used waxing tips

Item: 10155  
Item: 10255 (230V)

SPECIFICATIONS	
Voltage	AC 100-240V 50/60Hz 20W
Dimensions	6" W x 4.8" D x 2.8" H 152mm x 122mm x 71mm
Temperature Range	50°C-200°C 122°F-392°F
Weight	2.9 lbs. (1.3 kg)
<small>NOTE: Minimum voltage supply should not exceed ± 10% of the nominal value</small>	



### Digital Wax Carving Pencil – Single

#### In addition to standard features:

- Switch between two temperatures at the touch of a button
- Space-saving small footprint
- Compartment in rear of unit for easy storage of carving tips
- Includes two wax tips

Item: 10157  
Item: 10257 (230V)

SPECIFICATIONS	
Voltage	AC 110/220V 50/60Hz
Dimensions	3" W x 4.7" D x 2.5" H 76mm x 119mm x 64mm
Temperature Range	50°C-200°C 122°F-392°F
Weight	1.2 lbs. (0.54 kg)
<small>NOTE: Minimum voltage supply should not exceed ± 10% of the nominal value</small>	

## Digital Wax Carving Pencil – Dual

*In addition to standard features:*

- Easy-to-use control buttons
- Two-line circuit, designed for simultaneous work
- Convenient external cradles
- Comes complete with six of the most commonly used waxing tips

Item: 10150  
Item: 10250 (230V)

SPECIFICATIONS	
Voltage	AC 100/240V 50/60Hz 20W
Dimensions	4.7" W x 6.3" D x 2.9" H 119mm x 160mm x 74mm
Temperature Range	50°C-200°C 122°F-392°F
Weight	2.8 lbs. (1.27 kg)

NOTE: Minimum voltage supply should not exceed ± 10% of the nominal value



## The Heat Zone

- Heats wax tools fast without flame
- Reduces carbon build up on instruments
- Energy efficient as power is used only at the moment of heating
- No gas lines required, making unit mobile and operational anywhere an electric outlet is available
- Indicator light glows green for ready mode and red for when in use
- Safety light blinks if instrument is left heating too long
- Replaceable inserts keeps unit clean of wax residue (5 included)

Item: 10147  
Item: 10247 (230V)

SPECIFICATIONS	
Voltage	AC 100/120V 50/60Hz AC 220/240V 50/60Hz
Dimensions	3.5" W x 8.3" D x 3.5" H 89mm x 211mm x 89mm
Weight	1 lbs. (0.45 kg)

NOTE: Minimum voltage supply should not exceed ± 10% of the nominal value



## ACCESSORIES

### Tips (available individually)

- Item:
1. 10171 Offset Explorer
  2. 10172 PKT 1
  3. 10173 PKT 2
  4. 10174 Discoid
  5. 10175 Mini-Beavertail
  6. 10176 Beavertail
  7. 10177 Large Beavertail
  8. 10178 PKT 7
  9. 10179 1/2 Hollenback



### Electric Heating Tool

- Operates with a variety of milling machines and surveyors
  - Specially designed pencil top permits excellent heat conductivity
  - Five tips available for 0° (thin and thick gauge), 2°, 4° and 6°
- Item: 10159



### Replacement Pencil

- Item: 10163 Wax Carving Pencil  
Item: 10165 Touch Carving Pencil

# Precisely What You Need

## DIGITAL WAX POTS

### Standard Features

- Easy to operate
- Excellent insulation and heat conductivity
- Current and preset temperature readings are easily visible
- Heats wax quickly and maintains its liquid state
- Maintains desired temperature accurately and precisely
- 2-year warranty



### Digital Dipping Pot

#### In addition to standard features:

- Stores previous working temperature in memory, even if the unit has been turned off
- Easy selection of either Celsius or Fahrenheit temperature scales

Item: 10160

Item: 10260 (230V)

SPECIFICATIONS	
Voltage	AC 115/230V 50/60Hz
Dimensions	4.3" W x 5.4" D x 2.8" H 108mm x 137mm x 71mm
Temperature Range	50°C-120°C 122°F-248°F
Weight	1.9 lbs. (0.87 kg)

NOTE: Minimum voltage supply should not exceed  $\pm 10\%$  of the nominal value

### Wax Pot – Multi Well

#### In addition to standard features:

- Multi-well pot is capable of heating four different waxes simultaneously
- Built-in thermostat to accurately maintain desired temperature
- LED indicates when setpoint has been reached

Item: 10158

Item: 10258 (230V)

SPECIFICATIONS	
Voltage	AC 115/230V 50/60Hz
Dimensions	7.6" W x 4.5" D x 2.8" H 193mm x 114mm x 71mm
Temperature Range	50°C-120°C 122°F-248°F
Weight	2.5 lbs. (1.13 kg)

NOTE: Minimum voltage supply should not exceed  $\pm 10\%$  of the nominal value



more reasons. one source

Whip Mix Corporation • 361 Farmington Avenue • Louisville, KY 40209  
800-626-5651 • 502-637-1451 • Fax 502-634-4512 • www.whipmix.com



05/2015