




Compatible Printers with
VeriRESIN™ 3D Print Material



COMPATIBLE PRINTERS WITH WHIP MIX VeriRESIN™ 3D PRINT MATERIAL

PRINTERS	 VeriMODEL™			 VeriMODEL™ Ortho	 VeriTRAY™	 VeriGUIDE™	 VeriSPLINT™	 VeriCAST™	 VeriGUM™	 VeriDENT™ Try-In	 VeriFIT™ Splint
	IVORY	GOLDEN BROWN	GREY	WHITE	TURQUOISE	CLEAR	CLEAR	RED	PINK	A2	CLEAR
 VERIEKO (405nm)	YES	YES	YES	YES	YES	YES		YES	YES	YES	
 VERIBUILD (405nm)	YES	YES	YES	YES		YES		YES			
 ASIGA MAX (385/405nm)	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES
 ASIGA PRO 275 (385nm)	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES
 ASIGA PRO 4K (385nm)	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES
 ASIGA ULTRA (385nm)	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES



COMPATIBLE PRINTERS WITH WHIP MIX VeriRESIN™ 3D PRINT MATERIAL

PRINTERS	 VeriMODEL™			 VeriMODEL™ Ortho	 VeriTRAY™	 VeriGUIDE™	 VeriSPLINT™	 VeriCAST™	 VeriGUM™	 VeriDENT™ Try-In	 VeriFIT™ Splint
	IVORY	GOLDEN BROWN	GREY	WHITE	TURQUOISE	CLEAR	CLEAR	RED	PINK	A2	CLEAR
 ASIGA MAX 2 (385nm)	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES
 ACKURETTA SOL PLUS	YES	YES	YES		YES						
 ACKURETTA SOL (405nm)	YES	YES	YES	YES	YES	YES		YES	YES	YES	
 ACKURETTA FREESHape 120 (405nm)	YES	YES	YES	YES		YES		YES			
 ACKURETTA DIPLO (405nm)		YES		YES		YES		YES			
 RAPID SHAPE PRO 20 / PRO 30	YES	YES	YES	YES	YES	YES			YES	YES	YES



COMPATIBLE PRINTERS WITH WHIP MIX VeriRESIN™ 3D PRINT MATERIAL

PRINTERS	 VeriMODEL™			 VeriMODEL™ Ortho	 VeriTRAY™	 VeriGUIDE™	 VeriSPLINT™	 VeriCAST™	 VeriGUM™	 VeriDENT™ Try-In	 VeriFIT™ Splint
	IVORY	GOLDEN BROWN	GREY	WHITE	TURQUOISE	CLEAR	CLEAR	RED	PINK	A2	CLEAR
 RAPID SHAPE D50+	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES
 STRAUMANN P20+/P30+	YES	YES	YES		YES	YES			YES	YES	YES
 FLASHFORGE HUNTER (405nm)	YES	YES	YES	YES		YES	YES				
 MICROLAY VERSUS (385/405nm)		YES		YES		YES		YES			
 NEW PRO 3D NP1 (385/405nm)		YES		YES				YES			



COMPATIBLE PRINTERS WITH WHIP MIX VeriRESIN™ 3D PRINT MATERIAL

PRINTERS	 VeriMODEL™			 VeriMODEL™ Ortho	 VeriTRAY™	 VeriGUIDE™	 VeriSPLINT™	 VeriCAST™	 VeriGUM™	 VeriDENT™ Try-In	 VeriFIT™ Splint
	IVORY	GOLDEN BROWN	GREY	WHITE	TURQUOISE	CLEAR	CLEAR	RED	PINK	A2	CLEAR
 SHINING3D ACCUFAB-L4D DENTAL (405nm)	YES	YES	YES	YES	YES	YES		YES	YES	YES	YES
 PHROZEN SONIC 4K (405nm)	YES	YES	YES	YES	YES				YES		YES
 PHROZEN SONIC MINI (405nm)					YES				YES		
 PHROZEN SONIC XL 4K (405nm)	YES	YES	YES	YES	YES	YES			YES	YES	YES
 RAYSHAPE SHAPE 1+ DENTAL (405nm)	YES	YES	YES	YES	YES	YES			YES	YES	



COMPATIBLE PRINTERS WITH WHIP MIX VeriRESIN™ 3D PRINT MATERIAL

PRINTERS	 VeriMODEL™			 VeriMODEL™ Ortho	 VeriTRAY™	 VeriGUIDE™	 VeriSPLINT™	 VeriCAST™	 VeriGUM™	 VeriDENT™ Try-In	 VeriFIT™ Splint
	IVORY	GOLDEN BROWN	GREY	WHITE	TURQUOISE	CLEAR	CLEAR	RED	PINK	A2	CLEAR
 RAYSHAPE EDGE E1 (405nm)	YES	YES	YES	YES	YES	YES			YES	YES	
 UNIZ NBEE (405nm)	YES	YES	YES	YES					YES		
 MAKEX M- DENTAL		YES				YES		YES			
 ELEGOO SATURN 3 ULTRA (405nm)	YES	YES	YES			YES			YES		YES
 ELEGOO MARS 4 ULTRA (405nm)	YES	YES	YES			YES			YES		YES
 NEXA3D XiP						YES				YES	