YOUR SINGLE STOP FOR YOUR ORTHO MODELING NEEDS.

Elevate your dental experience with the Whip Mix® 3D-Printed Ortho and Presentation models.

- **Beautiful pure white color:** Color blends beautifully for any presentation model, even after post-processing
- **Ultimate accuracy coupled with smooth detailed surface:** Ensure all thermoformed devices fit perfectly.
- **High durability:** Robust enough to be used with a laser welder and handles carving without chipping.
- **Easy and fast post processing:** Experience efficiency with high production of thermoformed devices.
- **High quality surface:** Excellent for use with silicon-based separators and easy removal from thermoformed aligners.



"The white resin looks nice! I really like the color. Great for diagnostic models. Thanks."
DJ Zanon,CDT, Nu-Art Dental Inc

ORDER Information

Call Whip Mix Customer Service or your dealer to order.

800-626-5651

3 YEAR Shelf Life

385 405 nm

To learn more about our products, please scan the QR code.

Proudly formulated and manufactured in the United States.

VeriMODEL™ Ortho

385/405nm

Precision & Speed

Physical Properties

Physical Properties

Property	Value
Flexural Strength (ASTM D790)	90 MPa
Flexural Modulus (ASTM D790)	2200 MPa
Elongation at Break (ASTM D638)	14%
Viscosity	700 ± 100 cP at 24± 1 °C
Density	1.11g/ml



VeriMODEL™ Ortho White

A versatile material that is suitable for a diverse range of applications, including crafting precision orthodontic and presentation models. It can be used with silicone-based separators and vacuum formers and is robust enough for laser welding.

"What my dentists like about VeriMODEL White resin is that it is very accurate, detailed, and smooth, and it has a beautiful white color. They have complimented me about the models I print for them. I appreciate that I can rely on it for consistently beautiful and accurate orthodontic and presentation models"

Jane Hyatt, K&B Orthodontic Lab, Norwalk, Connecticut

INTRODUCING THE WHIP MIX® LINE OF 3D PRINT RESINS FOR THE DENTAL PROFESSIONAL

WHIP MIX CORPORATION, 361 Farmington Avenue, Louisville, KY USA 40209, www.whipmix.com



VeriDENT™ Try-In | VeriTRAY™ | VeriGUM™ | VeriGUIDE™ | VeriCAST™ | VeriMODEL™ Ivory | VeriMODEL™ Golden Brown | VeriMODEL™ Grey | VeriMODEL™ White | VeriFIT™ Splint